Certified quality

BioFix Bone Cements are CE marked through the European Directive 93/42/CEE, complying ISO 9001, ISO 13485 and ISO 5833 specifications for acrylic resin cements for surgical implants.

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LUBMEDICAL

BioFix Cement

BioFix Bone Cements are radiopaque acrylic compounds of spontaneous solidification intended to fix metal or plastic prosthesis to the live bone. Radiopacity is achieved by micro dispersed Barium Sulfate granules.

Types

There are 6 types of BioFix Bone Cements available: 3 to be manually applied and 3 for syringe application.

**MANUAL APPLICATION CEMENTS
(High viscosity)**

BioFix Manual Application Cements were developed to achieve:

- Effortless mixing: immediate homogeneity of dough.
- Uniform dough properties during a longer working period.
- Long working period with soft easy to use dough.
- Short final hardening stage, minimizing prosthesis micro movements during setting.
- Homogeneous radiopacity.
- Addition of antibiotic does not affect the physical, mixing, handling or insertion properties of the bone cement.
- Extended antibiotic release.
- Outstanding concentration at spongy and cortical bones.

Presentation:

The cements are available in blisters, containing one pouch with 40 g (ETO sterilised) powder , and one ampoule with 20 ml (ultra filtration sterilized) liquid.

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**SYRINGE APPLICATION CEMENTS
(Low viscosity)**

BioFix Syringe Application Bone Cements were developed to be used with any approved available syringe / gun system:

- Effortless mixing, providing immediate homogeneity of mixture.
- Easy application due to an extended low viscosity period.
- Reliable prosthesis placement and long lasting fixation.
- Short final hardening stage, minimizing prosthesis micro movements during setting.
- Homogeneous radiopacity.
- Addition of antibiotic does not affect the physical, mixing, handling or insertion properties of the bone cement.
- Extended antibiotic release.
- Outstanding concentration at spongy and cortical bones.

Presentation:

The cements are available in blisters, containing one pouch with 60 g (ETO sterilized) powder , and one ampoule with 30 ml (ultra filtration sterilized) liquid.

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BioFix PR Cement**Radiopaque surgical cement with Chlorophyll**

BioFix PR cements are used for bone cementations and fixations. Its use is advised in Orthopaedics and other surgical specialities, when fixation of prosthesis elements with sterile setting material is sought.

This line includes chlorophyll providing the cement its green colour which is easily distinguishable from the bone natural colour. The manual application materials demand a longer period of work and application and, consequently, a longer setting time.

Presentations

BioFix PR cement is used for manual application and is supplied in a package including a pouch with 40g of radiopaque powder sterilized by ethylene oxide and a blister with 20 ml of liquid sterilized by ultrafiltration.

BioFix PRG antibiotic loaded cement is used for manual application and is supplied in a package including a pouch with 40g of radiopaque powder with 0.50 g of gentamycin sterilized by ethylene oxide and a blister with 20 ml of liquid sterilized by ultrafiltration.

BioFix PR3 cement is used for gun application and is supplied in a package including a pouch with 60 g of radiopaque powder sterilized by ethylene oxide and a blister with 30 ml of liquid sterilized by ultrafiltration.

BioFix PR3G antibiotic loaded cement is used for gun application is supplied in a package including a pouch with 60 g of radiopaque powder with 0.75 g of gentamycin sterilized by ethylene oxide and a blister with 30 ml of liquid sterilized by ultrafiltration.

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BioFix Cranioplasty Cement**Surgical cement sterile for cranioplasty.****Specially formulated for:**

- Neurosurgery .
- Plastic surgery.
- Maxillo-facial surgery.
- ENT surgery.

BioFix Cranioplasty cement is intended for the repair of traumatic, tumoral or surgical bone defects.

Main advantages in comparison to orthopedic cements:

- Smaller PMMA particles.
- Improved manipulation and modeling time.
- Better impression taking.
- Smoother finishing.

Fixation methods:

- Non reabsorbable surgical sutures.
- Titanium wires or cables.
- Titanium plates / screws.

Presentation:

BioFix Cranioplasty is supplied in a double dose pack inside a blister sterilized by ethylene oxide.

Each dose consists of:

- 1 ampoule with 17 ml of liquid sterilized by ultrafiltration in a blister sterilized by ethylene oxide.
- 1 envelope with 30 g of acrylic powder sterilized by ethylene oxide

Each blister provides a triple sterility barrier.

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BioFix Cranioplasty - REF 880125

BioFix VTP Cement

Surgical cement for vertebroplasty

BioFix VTP, low viscosity acrylic bone cement has been specially designed to be used in percutaneous vertebroplasty, resulting in a compound with adequate working times combined with a high radiopacity to x-Rays.

Special features:

- Easy preparation and delivery .
- Low viscosity.
- No need of dyes, retardants, or solvents.
- Constant working times.
- Controlled exothermic reaction and outstanding mechanical properties.
- Enhanced view under x-Rays.
- Less exposure to x-Rays of Patients and OR Personnel.

Presentation:

BioFix VTP is available in a unit containing a pouch with 25 g radiopaque powder sterilized using ethylene oxide and an ampoule with 10 ml of liquid sterilized by ultrafiltration.

Each blister provides a triple sterility barrier .

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BioFix VTP Cement - REF 880825

BioFix VTP System

Cement application system for vertebroplasty

The BioFix VTP System is intended for the application of any Specific Vertebroplasty Cement. This set was designed only to be used with the "BioFix VTP GUN " since it will not fit to other application devices.

The BioFix VTP System includes:

- Vertebral body puncture needle with mandrin.
- Luer Lock m/f connector .
- 10 ml syringe.
- Syringe piston.
- Spatula.
- Mixing bowl.
- Surgical hammer.
- Suction cannula.

Presentation:

The BioFix VTP System case includes the instructions for use and a blister sterilized using ethylene oxide.

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
BioFix VTP Integral System

Integral cement application system for vertebroplasty cement.

The BioFix VTP Integral System is developed to apply any type of cement specifically designed for vertebroplasties. This set does not require any supplementary device since it includes a disposable application gun with a mounted syringe, which also simplifies the assembling of the set.

The BioFix VTP Integral System includes:

- 1 Vertebral body puncture cannula with mandrel
- 1 Luer Lock Male - Female connecting tube
- 1 10 ml syringe - gun assembly
- 1 Mixing bowl
- 1 Spatula
- 1 Cement suction cannula

Presentation:

The BioFix VTP Integral System is available in a case with the Instructions for Use and a double blister with the set sterilized with ethylene oxide.

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O 3.0 mm x 150 mm - REF 884029

BioMix VTP System

Vacuum mixer system for vertebroplasty cement

BioMix VTP is intended to mix a wide range of vertebroplasty cements.

The Vacuum Mixing device BioMix VTP for Vertebroplasty Cement includes:

- 1 Mixing Body
- 1 Body Cap with mixing rod
- 1 plastic Tube, 3 mm inner diameter and 2 m long
- 1 vapor-absorbing cartridge
- 1 Funnel
- 1 Raccord coupling for tubes (3 mm to 6 mm diameter)

Special features

- Cement can be prepared with the mixer directly connected to the Application Syringe. Prevents exposure to fumes.
- Reduces entrapment of air bubbles during mixing.
- Avoids transfer of cement from mixer into syringe as vertebroplasty cement is directly aspirated into the Delivery Syringe.
- Compatible to most vertebroplasty systems.

Presentation

BioMix VTP Vacuum device is available in a case containing the instructions for use and a blister with BioMix VTP elements within another blister, thus assuring a double protecting sterile barrier. Sterilization is performed with ethylene oxide.

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BioMix VTP - REF 884016

BioMix V System**Vacuum mixer / delivery system**

Vacuum mixer and application system for bone cement.

Mixing and application device for orthopedic cement
BioMix V includes:

- 1 Mixing Body
- 1 Body Cap with mixing rod
- 1 plastic Tube, 3 mm inner diameter and 2 m long
- 1 vapor-absorbing cartridge
- 1 Funnel
- 1 Raccord coupling for tubes (3 mm to 6 mm diameter)
- 1 Pressurizer

Special features

- Reduces entrapment of air bubbles.
- Improves mechanical properties of cements.
- Trouble-free cement application.
- Prevents exposure to fumes.
- Cannula is pre-cut to fit cement delivery in knee and/or hip surgeries.
- Fits current application guns.
- Central application cannula allows cementation in narrow or confined surgical fields.

Presentation

BioMix V system is available in a case containing the instructions for use and a blister with BioMix V elements within another blister, thus assuring a double protecting sterile barrier. Sterilization is performed with ethylene oxide.

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- 1 Cartridge with plunger
- 1 Nozzle adapter
- 1 Long Nozzle
- 1 Short Nozzle
- 1 Mixing Bowl
- 1 Spatula
- 1 Funnel
- 1 Pressurizer (optional)

Special features

- Once prepared, cement is poured into application cartridge.
- Fully disposable.
- Two lengths of cannulas: knee and hip cementation.
- To be used with BioFix GUN device.
- Cost-efficient.
- Easy operation.

Presentation

BioFix 3 Syringe is available in a box with the Instructions for Use and a double blister with the bone cement dispenser system sterilised by ethylene oxide gas.

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BioFix spacer - Hip

BioFix Hip Spacer is intended to be temporarily implanted as an aid in the treatment of prosthetic hip infections, by means of a two stages surgical technique.

Characteristics:

5 models built with an AISI 316L surgical grade stainless steel core.
High resistance bone cement body, antibiotic loaded (Gentamicin).

Advantages:

- Efficient local (Gentamicin) release.
- Metallic core assures resistance to physiological and mechanical stress.
- Prevents proximal femoral migration, facilitating the surgical approach during the final revision.
- Preserves hip muscular function.
- Simplifies the surgical act.
- Early patient rehabilitation due to premature mobilization.
- Improved quality of life between surgical steps.

	Stem Lengths (mm)	Stem Type	Head Diameters (mm)	Angle°	REF
40C	127	Charnley	40	133	881150
48C	130		48		881151
48M	120	Müller	48		881153
56C	130	Charnley	56		881152
56M	127	Müller			881154

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**BioFix Spacer - Knee**

BioFix Knee Spacer is intended to be temporarily implanted as an aid in the treatment of prosthetic knee infections, by means of a Two Stages surgical technique.

Characteristics:

3 Models, completely built in thermocurable acrylic cement, Antibiotic (Gentamicin) loaded.

Advantages:

- Efficient local (Gentamicin) release.
- Maintains joint space.
- Preserves joint movement. Simplifies surgical approach, since patellar tendon is free from adhesions produced by immobilization.
- Maintains extensor mechanism and preserves ligament functionality due to use.
- Preserves Lateral Ligaments functional.
- Better mobilization range than with block spacers.
- Reduces functional recovery period after a definitive revision surgery.
- Improved quality of life between surgical steps.

	Articular Surface Diameter (mm)	REF
58	54	881155
65		881156
79	68	under design

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Test Prostheses BioFix**Set of metallic test prostheses**

During surgery, before opening the spacer, the selected size should be confirmed by placing a test prosthesis.

Presentation

- **Hip test prostheses:** perforated aluminium container for sterilization, with stainless steel test prostheses corresponding to the five types of hip prosthesis. All of them are held in place with silicone brackets.
- **Knee test prostheses:** perforated aluminium container for sterilization with stainless steel tibial and femoral test prostheses corresponding to the two types of knee prosthesis. All of them are held in place with silicone brackets.

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**BioBowl****Mixing kit**

The BioBowl includes: 1 Mixing Bowl / 1 Spatula

Presentation

BioBowl is available in a double blister with the mixing kit sterilised by ethylene oxide gas.

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**BioFix VTP Cannula****Vertebral body puncture cannula with handle and mandrel**

This cannula is used for puncturing the vertebral body during Vertebroplasty.

Presentation:

Each container includes a blister with a cannula sterilized with ethylene oxide.

Ř 2,4 mm x 100 mm - REF 884024
Ř 3,0 mm x 100 mm - REF 884025
Ř 3,0 mm x 150 mm - REF 884026

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**BioMix V GUN****Reusable injecting device to inject orthopaedic cement to the femoral canal with Mixing and application Set BioMix V.****Characteristics and advantages:**

- Easy operation.
 - Ergonomic design.
 - Autoclave, dry heat and ETO sterilization.
 - Immediate interruption of cement injection in case of undesirable events.
 - Made of anodised aluminium or stainless steel.
- BioMix V Gun is supplied in a sterilizable box.

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BioFix VTP GUN**Reusable injecting device for vertebroplasty cement.**

Made of anodised aluminium and stainless steel.

BioFix VTP Gun is supplied in a sterilizable box and has been specially designed to be use with BioFix VTP System in percutaneous vertebroplasties.

Characteristics and advantages:

- Easy operation.
- Ergonomic design.
- Autoclave, dry heat or ETO sterilization.
- Exact cement injection under low and high pressure.
- Immediate interruption of cement injection in case of undesirable events.

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
**BioFix GUN****Reusable device for the injection of orthopaedic cement, specially designed to be use with BioFix 3 Syringe.****Characteristics and advantages:**

- Very easy to clean.
 - Totally made of stainless steel.
 - No wearing out after successive sterilizations of any type.
- BioFix GUN is supplied in a sterilizable box.

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