

Videodrucker | Imprimante vidéo couleur
 Video printer | Impresora video color
 Stampante a colori



Sony UP-25MD **tk 753-0170**

Videodrucker der neuesten Generation
Video printer of the latest generation

Technische Daten

• TV-System	PAL / NTSC
• Auflösung	423 dpi
• Farben	> 16,7 Mio.
• HDTV IN/OUT - Signale	1080i / 720p
• Eingänge	RGB, S-Video, Y/C, BNC
• Stromversorgung	100 - 240 V AC, 50/60 Hz
• Stromaufnahme	1,7A - 1,0A
• Betriebstemperatur	+5°C bis +35°C
• Luftfeuchtigkeit	20% - 80%
• Abmessungen (B x T x H)	212 x 398 x 98 mm
• Gewicht	5,7 kg
• Papierformat	A6-S (100 x 90 mm) A6-L (144 x 100 mm)

Im Lieferumfang enthalten:

- CD-ROM (mit Bedienungsanleitung)
- Papierfach
- Reinigungspatrone

Zubehör
Accessories

tk 753-0178

UPC-21S Druckpaket für 240 Bilder (100 x 90 mm)
UPC-21S Printpackage for 240 images (100 x 90 mm)

Technical data

• TV-System	PAL / NTSC
• Resolution	423 dpi
• Colors	> 16.7 Mio.
• HDTV IN/OUT - Signals	1080i / 720p
• Inputs	RGB, S-Video, Y/C, BNC
• Power supply	100 - 240 V AC, 50/60 Hz
• Power consumption	1.7A - 1.0A
• Operating temperature	+5°C to +35°C
• Humidity	20% - 80%
• Dimensions (W x D x H)	212 x 398 x 98 mm
• Weight	5.7 kg
• Paper size	A6-S (100 x 90 mm) A6-L (144 x 100 mm)

Included in delivery

- CD-ROM (with operating manual)
- Paper tray
- Cleaning cartridge

tk 753-0179

UPC-21L Druckpaket für 200 Bilder (144 x 100 mm)
UPC-21L Printpackage for 200 images (144 x 100 mm)